

Oosafe®

PLASTIC WARE FOR IVD USE



CE Manufactured in accordance
with EN ISO 13485:2016

Oosafe® PLASTIC WARE FOR IVD USE



TOP QUALITY PLASTICWARE USED BY SCIENTISTS ALL OVER THE WORLD

- Designed for wide range of applications
- Tested for each LOT:
 - One cell stage Mouse Embryo Assay (MEA > 80% Blastocyst)
 - Non-pyrogenic, Endotoxin (LAL < 0.03 EU/device)
 - Human Sperm Survival and Motility (SMI ≥ 0.75)
- Surface treated options for a consistent hydrophilic surface
- Packaged in peel or tear-open packaging for user convenience
- Easy Grip and Label Area Design
- Sterility assurance level (Irradiation) 10⁻⁶
- CE marked according to 98/79/EC Directive on *In Vitro* Diagnostic Medical Devices
- Made of USP Class VI tested Raw Materials
 - Polystyrene (PS)
 - Polypropylene (PP)
 - Polyethylene (PE)
- Vacuum packed to prevent scratches



Oosafe® 4 Well Dish TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-FW03

*Cut corner design ensures proper orientation
and prevents mix-ups while handling*
Well volume 1.8 mL
Volume between wells 7 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 4 Well Dish NON-TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-FW04

*Cut corner design ensures proper orientation
and prevents mix-ups while handling*
Well volume 1.8 mL
Volume between wells 7 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 6 Well Dish with Straw Holder

SINGLEPACK PS

1 pc/pack
240 pcs/case



Product reference:
OOPW-SW02

*Conical flat bottom. Hands-free design: the tip of the straw
and the bottom of the well are at the same focus level*
Well volume 0.4 mL
Volume between wells 3 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 6 Well Dish with Straw Holder

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-SW03

*Conical flat bottom. Hands-free design: the tip of the straw
and the bottom of the well are at the same focus level*
Well volume 0.4 mL
Volume between wells 3 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Center Well Dish with Two Compartments



LABEL AREA GRIP
SINGLEPACK PS

1 pc/pack
240 pcs/case



Product reference:
OOPW-CW04

Multi-purpose dish for cell handling applications
Center well volume 3 mL
Each compartment volume 5 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Center Well Dish with Two Compartments



LABEL AREA GRIP
MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-CW05

Multi-purpose dish for cell handling applications
Center well volume 3 mL
Each compartment volume 5 mL

Oosafe® Center Well Dish



LABEL AREA GRIP
MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-CW08

New design with one outer chamber
Center well volume 3 mL
Outer compartment volume 10 mL

Oosafe® 50 mm Dish



LABEL AREA GRIP
SINGLEPACK PS

1 pc/pack
240 pcs/case



Product reference:
OOPW-IC02

Extra thin observation area for better image quality
Thickness of observation area: 0.9 mm
Thickness of outer ring: 1 mm
Volume 12 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® 50 mm Dish



LABEL AREA GRIP
MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-IC03

Extra thin observation area for better image quality
Thickness of observation area: 0.9 mm
Thickness of outer ring: 1 mm
Volume 12 mL

Oosafe® 60 mm Dish,



LABEL AREA GRIP
MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-ST03

Multi-purpose dish for cell culture
and microscopy applications
Volume 25 mL

Oosafe® 35 mm Dish, High Wall



LABEL AREA GRIP
MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-TF03

Multi-purpose dish for cell culture and microscopy
applications, high wall with enough space
for labels between the grips
Volume 8 mL

Oosafe® 100 mm Dish



MULTIPACK PS

10 pcs/pack
250 pcs/case



Product reference:
OOPW-HD10

Multi-purpose dish for cell culture
and microscopy applications
Volume 105 mL

Oosafe® Synergy Dish

TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:

OOPW-SD04

Unique design with micro-wells for group culture

Volume of wash well 56 µl

Volume of group culture well 160 µl

Total volume 15 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Andrology Tube, 5 mL

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:

OOPW-AT10

Push-on cap with dual position to enable venting

Minimum volume measured 0.5 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® OPU Tube, 14 mL

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:

OOPW-OT10

Push-on cap with dual position to enable venting

Minimum volume measured 0.1 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Centrifuge Tube, 15 mL

SINGLEPACK PP

1 pc/pack
300 pcs/case



Product reference:

OOPW-CT01

Suited for applications requiring greater

thermal and chemical stability

Minimum volume measured 0.1 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® PS Centrifuge Tube, 15 mL

MULTIPACK PS

5 pcs/pack
500 pcs/case



Product reference:

OOPW-CT05

Suited for procedures requiring high optical clarity
Minimum volume measured 0.1 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® PP Centrifuge Tube, 50 mL

MULTIPACK PP

25 pcs/pack
500 pcs/case



Product reference:

OOPW-CT25

Ideal for larger volume centrifugation
Minimum samples volume measured 5 mL

Oosafe® Sperm Collection Cup, 80 mL

SINGLEPACK PP

1 pc/pack
70 pcs/case



Product reference:

OOPW-SC01

Wide opening for easy sampling
Minimum volume measured 2 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Flask, 50 mL

TREATED SURFACE

MULTIPACK PS

10 pcs/pack
160 pcs/case



Product reference:

OOPW-FL10

Bent neck allows to easily access all corners of the flask
Minimum volume measured 5 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® PIPETTE COLLECTION



Oosafe® Pasteur Pipette, 1 mL

SINGLEPACK PE

1 pc/pack
500 pcs/case



Product reference:
OOPW-PP01

For quick and safe transfer of media and samples

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Pasteur Pipette, 3 mL

SINGLEPACK PE

1 pc/pack
500 pcs/case



Product reference:
OOPW-PP03

For quick and safe transfer of media and samples

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® 1 mL in 1/100 mL Serological Pipette

SINGLEPACK PS

1 pc/pack
500 pcs/case



Product reference:
OOPW-SP01

Negative graduations for extra capacity

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® 2 mL in 1/100 mL Serological Pipette

SINGLEPACK PS

1 pc/pack
500 pcs/case



Product reference:
OOPW-SP02

Negative graduations for extra capacity

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® 5 mL in 1/10 mL Serological Pipette

SINGLEPACK PS

1 pc/pack
250 pcs/case



Product reference:
OOPW-SP05

*Negative graduations for extra capacity
and reverse graduations for user convenience*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® 10 mL in 1/10 mL Serological Pipette

SINGLEPACK PS

1 pc/pack
200 pcs/case



Product reference:
OOPW-SP10

*Negative graduations for extra capacity
and reverse graduations for user convenience*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

PRODUCTS AVAILABLE UPON REQUEST *



**Oosafe® Center Well Dish
with Two Compartments**
ALL AROUND GRIP

MULTIPACK **PS**

10 pcs/pack
500 pcs/case

*All around grip for convenient and
secure handling*

Product reference:
OOPW-CW03



Oosafe® 60 mm Dish,
LABEL AREA GRIP
TREATED SURFACE

MULTIPACK **PS**

10 pcs/pack
500 pcs/case

*Surface treated for a consistent
hydrophilic surface*

Product reference:
OOPW-ST05



**Oosafe® 35 mm Dish,
High Wall**
LABEL AREA GRIP
TREATED SURFACE

MULTIPACK **PS**

10 pcs/pack
500 pcs/case

*Surface treated for a consistent
hydrophilic surface*

Product reference:
OOPW-TF05

*Minimum purchase quantity of 1000 cases.



All plasticware products are packed in a secure white cardboard box.



PLASTIC WARE

FOR IVD USE

ORDER INFORMATION

Product reference	Product name	Surface treatment	Packaging
OOPW-FW03	Oosafe® 4 Well Dish	Treated	MULTIPACK
OOPW-FW04	Oosafe® 4 Well Dish	Non-treated	MULTIPACK
OOPW-SW02	Oosafe® 6 Well Dish with Straw Holder	Non-treated	SINGLEPACK
OOPW-SW03	Oosafe® 6 Well Dish with Straw Holder	Non-treated	MULTIPACK
OOPW-CW04	Oosafe® Center Well Dish with two Compartments, Label Area Grip	Non-treated	SINGLEPACK
OOPW-CW05	Oosafe® Center Well Dish with two Compartments, Label Area Grip	Non-treated	MULTIPACK
OOPW-CW08	Oosafe® Center Well Dish	Non-treated	MULTIPACK
OOPW-IC02	Oosafe® 50 mm Dish, Label Area Grip	Non-treated	SINGLEPACK
OOPW-IC03	Oosafe® 50 mm Dish, Label Area Grip	Non-treated	MULTIPACK
OOPW-ST03	Oosafe® 60mm Dish, Label Area Grip	Non-treated	MULTIPACK
OOPW-TF03	Oosafe® 35mm Dish, High Wall	Non-treated	MULTIPACK
OOPW-HD10	Oosafe® 100mm Dish	Non-treated	MULTIPACK
OOPW-SD04	Oosafe® Synergy Dish	Treated	MULTIPACK
OOPW-AT10	Oosafe® Andrology Tube, 5mL	Non-treated	MULTIPACK
OOPW-OT10	Oosafe® OPU Tube, 14 mL	Non-treated	MULTIPACK
OOPW-CT01	Oosafe® Centrifuge Tube, 15 mL	Non-treated	SINGLEPACK
OOPW-CT05	Oosafe® PS Centrifuge Tube, 15 mL	Non-treated	MULTIPACK
OOPW-CT25	Oosafe® PP Centrifuge Tube, 50 mL	Non-treated	MULTIPACK
OOPW-SC01	Oosafe® Sperm Collection Cup, 80 mL	Non-treated	SINGLEPACK
OOPW-FL10	Oosafe® Flask, 50 mL	Treated	MULTIPACK
OOPW-PP01	Oosafe® Pasteur Pipette, 1 mL	Non-treated	SINGLEPACK
OOPW-PP03	Oosafe® Pasteur Pipette, 3 mL	Non-treated	SINGLEPACK
OOPW-SP01	Oosafe® 1 mL in 1/100 mL Serological Pipette	Non-treated	SINGLEPACK
OOPW-SP02	Oosafe® 2 mL in 1/100 mL Serological Pipette	Non-treated	SINGLEPACK
OOPW-SP05	Oosafe® 5 mL in 1/10 mL Serological Pipette	Non-treated	SINGLEPACK
OOPW-SP10	Oosafe® 10 mL in 1/10 mL Serological Pipette	Non-treated	SINGLEPACK
Products available upon request. Please contact our Sales department: sales@sparmed.dk			
OOPW-CW03	Oosafe® Center Well Dish with two Compartments	Non-treated	MULTIPACK
OOPW-ST05	Oosafe® 60mm Dish, Label Area Grip	Treated	MULTIPACK
OOPW-TF05	Oosafe® 35mm Dish, High Wall	Treated	MULTIPACK

Oosafe[®]
**PLASTIC
WARE**
FOR IVD USE



SparMED ApS
Address: Ryttermarken 2,
3520 Farum, DENMARK
Phone number: +45 39-40-25-03



MK_4.3.B2 rev. 07.01.2021